

ThermoFisher NanoDrop 8000



To request a quote contact: McClernonLLC@gmail.com

Description of Item: Year 2014 NanoDrop 8000 Spectrophotometer system in excellent condition and comes with new ThermoFisher PM service August 2018.

Comes with:

- NanoDrop 8000 Spectrophotometer Year 2014
- HP ProBook 4535s laptop with WIN 7 OS loaded with newest v2.3.2 software
- Comm USB Cable
- New ThermoFisher PM Service August 2018
- 6 Month Warranty

Perfect for high-throughput amplicon analysis!

- Improved productivity, capable of analyzing one to eight 1 μ L samples at one time
- Measure up to 96 samples in less than six minutes
- Innovative, user-friendly software
- Sample Position Illuminator reduces error by keeping track of the samples to be measured
- Full-spectral data displayed for each sample (220 – 750 nm)
- Measure nucleic acid concentration and purity (2.5 – 3,700 ng/ μ L for dsDNA)
- Measure protein concentration and purity (0.15 – 100 mg/mL for BSA)

[CLICK HERE TO VIEW ADDITIONAL ITEMS](#)